

# GroMed™

## 3 grams Soybean - Casein Digest Medium

### Catalog #GM3000

For the cultivation of microorganisms. Not for patient use.

Non-sterile.

Store at 2 to 30° C.

1. Combine contents of packet (approximately 3gm) with 100mL of non-bacteriostatic water to make a 3% solution.
2. Heat and stir as needed to dissolve completely.
3. Yields approximately 100mL of non-sterile Soybean-Casein Digest Broth (Tryptic Soy Broth).
4. Liquid media can be inoculated with specific microorganisms or set out over night in an uncovered wide-mouth container.
5. Optional: Use in conjunction with PATT™ kits, #GM7020 or #GM7030, to perform a high-risk media fill test.
  - a. Follow all the steps described in the PATT or PATT 2 directions for use.
  - b. Dissolve 3 grams of Soybean-Casein Digest Media, GroMed #GM3000, in 100mL of non-sterile water.
  - c. Draw 10 mL of the non-sterile TSB into a 20 mL syringe. Attach a 0.2µ syringe filter and sterile needle to the syringe.
  - d. Inject 10 mL of sterilized TSB into the bag. Label and incubate as described in the PATT or PATT 2 directions for use.

**IRRITANT:** Dry Powder may be irritating to eyes, skin, and respiratory system.  
Avoid contact or inhalation.

## **Q.I.medical, inc**

Grass Valley, CA 95945  
(530) 272-8700 • FAX (530) 272-8702  
www.qimedical.com

Q.I. Medical, Inc.  
Price List



DFU020 E 06/14

