


Q.I.medical, inc



Price List

Effective: January 1, 2005

Catalog #	Description	Case Contents	Cost/ea.	Cost/case
-----------	-------------	---------------	----------	-----------

EnviroTest™ Test Kits

ET1000	Double sided growth media paddles for surface, air, or fingertip contamination testing, formulated to inactivate many antimicrobial compounds.	10 Testers	\$3.90	\$39.00
--------	--	------------	--------	---------

Gro Med™

Aseptic Technique Training and Testing Aids

NEW	GM0100	TSB Growth media, 10ml tube screw cap and needle access port.	20 Tubes	\$2.45	\$49.00
	GM0200	TSB Growth media, 20ml vial, needle/spike septum.	20 Vials	\$3.15	\$63.00
	GM1400	TSB Growth media, 100ml bag with needle access port and spike port.	10 Bags	\$5.60	\$56.00
	GM1500	TSB Growth media, 100ml bag with Male Luer Lock port and spike port.	10 Bags	\$6.25	\$62.50
	GM2000	TSB Growth media, 100ml vial with septum for needle or dispensing pin access.	10 Vials	\$4.79	\$47.90
NEW	GM3000	Soybean-Casein Digest Media, approximately 3 grams of nonsterile powder in foil pouch, makes 100mL of TSB media	20 Pouches	\$3.60	\$72.00
NEW	GM4000	PDA, Potato Dextrose Agar, 10ml tubes, test for fungal contamination	20 Tubes	\$3.25	\$65.00
	GM5000	TSB Growth media, 500ml bag, 1X strength, needle access port and spike port.	10 Bags	\$9.75	\$97.50
	GM6005	5cc TSB Growth media in 12cc luer lock syringe.	10 Syringes	\$2.45	\$24.50
NEW	GM6105	5cc FTG (Fluid Thioglycollate) growth media in 12cc luer lock syringe.	10 Syringes	\$2.45	\$24.50
	GM7020	PATT, Personal Aseptic Technique test kit, 100ml partially filled minibags, test procedure result log.	5 Vials 5 Bags	\$9.99	\$49.95
	GM7030	PATT2, Personal Aseptic Technique test kit, 3ml ampules, 20ml vials, 100ml partially filled minibags, test procedure and results log.	5 Ampules 5 Vials 5 Bags	\$13.95	\$69.75

Quick Test™ System

Full filtration for all LVPs, including emulsions

QT1000	Contamination tester for use with: • ADComp™ Compounder • Automix® Compounder (Baxter) • Autocomp 6™ Compounder • Exacta-Mix™ Compounder (Baxa) • MicroMacro™ Compounder (Baxa) • Nutrimix™ Compounder (Abbott)	10 Testers 10 100ml bags of media	\$19.50	\$195.00
QT1100	Contamination tester with universal spike and 18" extension set.	10 Testers 10 100ml bags of media 10 16 ga. needles	\$22.00	\$220.00
QT2000	Contamination tester for use with: • Any pump set with male luer connector • Repeater™ Pump (Baxa) • Hyperformer™ Compounder (Braun) • PharmAssist™ Pump (Excelsior)	10 Tester 10 100ml bags of media	\$20.15	\$201.50

Q.I.medical, inc

P.O. Box 1329, NevadaCity, CA 95959
530-265-4820 Fax 530-265-9416
qimed@qimedical.com
www.qimedical.com

Price List *continued*

Catalog #	Description	Case Contents	Cost/ea.	Cost/case
Tuff Test™ System				
Aliquot test for emulsions and suspensions				
NEW GM8000	Tuff Test™ contamination tester, 70ml TSB and 2 agar paddles in sterile plastic bottle, needle access port.	4 Testers	\$13.50	\$54.00
QT Junior™ System				
Full filtration for all LVPs, excluding emulsions				
TJ3100	Contamination tester with .22µ membrane universal spike inlet and male luer outlet	20 Testers 20 20ml vials of media	\$10.65	\$213.00
TJ3200	Contamination tester with .22µ membrane universal spike inlet and spike septum outlet	20 Testers 20 20ml bags of media	\$10.65	\$213.00
QT Micro™ Systems				
Full filtration for small volume sterile solution				
TM6000	0.22µ Contamination tester, female luer lock inlet female luer lock outlet.	10 Testers 10 5cc syringes of media	\$7.85	\$78.50
TM6100	0.22µ Contamination tester, female luer lock inlet, male luer slip outlet.	10 Tester 10 5cc syringes of media	\$7.85	\$78.50
PyroTest™ Test Kits				
Endotoxin Detection Test Kits				
PT5000	Endotoxin test for compounded sterile drug products.	5 Tests	\$30.70	\$153.50
INC004	PyroTest™ vial block, holds 4 vials upright in incubator, aluminum (use with QL 10-140 or QL12-140 incubators).	1 Vial Block	\$16.15	N/A
GroMed™				
Incubators				
QL10-140	Incubator 13"W x 11"D x 15"H, 3 shelves, dial thermometer.	1 Incubator	\$399.00	N/A
QL12-140	Large incubator 19"W x 13"D x 21"H, 3 shelves, dial thermometer.	1 Incubator	\$499.00	N/A
Wizard™ Syringe Filters				
DF010H	0.10µ filter, 25mm, sterile, female luer lock inlet, male luer lock outlet.	50 Filters per dispenser box, 4 boxes per case	\$1.66	Box \$83.00 Case \$332.00
DF022H	0.22µ filter, 25mm, sterile, female luer lock inlet, male luer lock outlet.	50 Filters per dispenser box, 4 boxes per case	\$1.55	Box \$77.50 Case \$310.00
NEW PG005	Bubble Point Tester, 0-100 psi range	1 Tester	\$73.60	N/A
Safe-Fill™ Vial Adaptors				
SF100	Universal vial adaptor, female luer slip outlet, no locking flanges.	50 Adaptors per box	\$0.85	\$42.50
SF200S	Small vial adaptor, female luer slip outlet, locking flanges.	50 Adaptors per box	\$0.90	\$45.00
SF200L	Large vial adaptor, female luer slip outlet locking flanges.	50 Adaptors per box	\$0.90	\$45.00

Q.I.medical, inc

P.O. Box 1329, Nevada City, CA 95959
 530-265-4820 Fax 530-265-9416
 qimed@qimedical.com
 www.qimedical.com