



**Price List**

**Effective: January 1, 2017**

Catalog #	Description	Case Contents	Cost/ea.	Cost/case
-----------	-------------	---------------	----------	-----------

**EnviroTest® Test Kits**

ET1000	Double sided TSA (trypticase soy agar) growth media paddles for surface or glove fingertip sampling, formulated to inactivate many antimicrobial compounds.	10 Testers	\$4.85	\$48.50
ET3000	Double sided MEA (malt extract agar) growth media paddles for surface or glove fingertip sampling, with added antibiotics to prevent bacterial overgrowth.	10 Testers	\$4.85	\$48.50

**GroMed®**

**Aseptic Technique Training and Testing Aids**

GM0030	TSB Growth media, 3mL ampule	20 Ampules	\$4.55	\$91.00
GM0100	TSB Growth media, 10mL tube screw cap and needle access port.	20 Tubes	\$3.80	\$76.00
GM0101	FTM Growth media, 10mL tube screw cap and needle access port.	20 Tubes	\$3.80	\$76.00
GM0200	TSB Growth media, 20mL vial, needle septum.	20 Vials	\$4.55	\$91.00
GM0201	FTM Growth media, 20mL vial, needle/spike septum	10 Vials	\$4.55	\$45.50
GM0203	Fluid D, 20mL vial, needle/spike septum	20 Vials	\$4.55	\$91.00
GM0252	TSB Growth media, 25mL vial, 2X, needle septum	10 Vials	\$6.09	\$60.90
GM1400	TSB Growth media, 100mL bag with covered needle port and spike port.	10 Bags	\$8.49	\$84.90
GM1500	TSB Growth media, 100mL bag with Male Luer Lock port and spike port.	10 Bags	\$9.59	\$95.90
GM2000	TSB Growth media, 100mL vial with septum for needle or dispensing pin access.	10 Vials	\$7.29	\$72.90
GM2001	FTM Growth media, 100mL vial with septum for needle or dispensing pin access.	10 Vials	\$7.19	\$71.90
GM3000	Soybean-Casein Digest Media, approximately 3 grams of nonsterile powder in foil pouch, makes 100mL of TSB media	20 Pouches	\$4.39	\$87.80
GM4000	PDA, Potato Dextrose Agar, 10mL tubes, test for fungal contamination	20 Tubes	\$4.19	\$83.80
GM5000	TSB Growth media, 500mL bag, 1X strength, covered needle port and spike port.	10 Bags	\$12.79	\$127.90
GM6005	5cc TSB Growth media in 12cc luer lock syringe.	10 Syringes	\$3.09	\$30.90
GM6105	5cc FTM growth media in 12cc luer lock syringe.	10 Syringes	\$5.19	\$51.90
GM7020	PATT, Personal Aseptic Technique test kit, 20mL vials, 100mL partially filled minibags, test procedure results log.	5 Vials 5 Bags	\$14.06	\$70.30
GM7030	PATT2, Personal Aseptic Technique test kit, 3mL ampules, 20mL vials, 100mL partially filled minibags, test procedure and results log.	5 Ampules 5 Bags, 5 Vials	\$19.80	\$99.00
EV0200	20mL sterile empty vial, not for patient use.	20 Vials	\$3.49	\$69.80
EV2000	100mL sterile empty vial, not for patient use.	10 Vials	\$5.09	\$50.92

**Quick Test® System**

**Full filtration for all LVPs, including emulsions**

QT1000	Contamination tester for use with most compounders	10 Testers/media	\$30.99	\$309.90
QT1100	Contamination tester with universal spike and 18" extension set.	10 Testers/media	\$33.99	\$339.90
QT2000	Contamination tester for use with any pump set with a male luer connector.	10 Testers/media	\$31.99	\$319.90
QT3000	Contamination tester for use with: Pinnacle compounder (BBraun)	10 Testers/media	\$30.99	\$309.90

Catalog #	Description	Case Contents	Cost/ea.	Cost/case
<b>TuffTest 2® System</b>				
<b>Aliquot test for emulsions and suspensions</b>				
GM8100	TuffTest 2™ contamination tester, 100mL TSB with TSA in sterile plastic bottle, needle access port.	4 Testers	\$19.04	\$76.18
GM8101	TuffTest 2™ contamination tester, 100mL FTM with Clostridial Agar in sterile plastic bottle, needle access port.	4 Testers	\$19.04	\$76.18
<b>QT Junior® System</b>				
<b>Full filtration for all LVPs, excluding emulsions</b>				
TJ3100	Contamination tester with .22µ membrane universal spike inlet and male luer outlet	20 Testers 20 20ml vials of media	\$15.95	\$319.09
TJ3200	Contamination tester with .22µ membrane universal spike inlet and spike septum outlet	20 Testers 20 20ml vials of media	\$15.95	\$319.09
<b>QT Micro® Systems</b>				
<b>Full filtration for small volume sterile solutions</b>				
TM6000	0.22µ Contamination tester, female luer lock inlet	10 Testers 10 5cc syringes of media	\$12.29	\$122.90
TM6001	0.22µ Contamination tester, FTM media		\$13.29	\$132.90
TM6100	0.22µ Contamination tester, female luer lock inlet, male luer slip outlet, TSB Media	10 Testers 10 5cc syringes of media	\$12.29	\$122.90
TM6101	0.22µ Contamination tester, FTM media		\$13.29	\$132.90
<b>ChemoTest®</b>				
<b>Hazardous Drug Handling Test Kits</b>				
CT9530	Validate the technique of health care workers who manipulate hazardous drugs. Uses fluorescein dye and a UV light to detect breaks in technique. 3 mL and 50 mL vials.	5 of each vial	\$29.99	\$149.95
CT9000	Mini Ultra Violet light, for use with ChemoTest #9530. Uses 4 AA batteries, not included.	1 Light	\$62.49	N/A
<b>GroMed®</b>				
<b>Incubators</b>				
QL10-140	Incubator 13"W x 11"D x 15"H, 3 shelves, dial thermometer.	1 Incubator	\$545.99	N/A
QL12-140	Large incubator 19"W x 13"D x 21"H, 3 shelves, dial thermometer.	1 Incubator	\$649.99	N/A
<b>Wizard™ Syringe Filters</b>				
HP1002	0.22µ filter, 25mm, sterile, female luer lock inlet, male luer lock outlet.	50 Filters per box	\$2.33	Box \$116.50
PG005	Bubble Point Tester, 0-100 psi range	1 Tester	\$199.99	N/A
<b>SafeFill® Vial Adaptors</b>				
SF100	Universal vial adaptor, female luer slip outlet, no locking flanges.	50 Adaptors per box	\$1.48	\$74.16
SF200S	Small vial adaptor, female luer slip outlet, locking flanges.	50 Adaptors per box	\$1.53	\$76.74
SF200L	Large vial adaptor, female luer slip outlet locking flanges.	50 Adaptors per box	\$1.53	\$76.74

